

MR 30 for SIGNA™

ARIA Monitoring Consensus Protocol

The following details define the GE HealthCare SIGNA™ MR scan protocol that meets the American Society for Neuroradiology (ASNR) Consensus criteria for the evaluation of Amyloid Related Imaging Abnormalities (ARIA)¹.


Essential Series

- 3 Plane Localizer
- 2D Axial T2 FLAIR
- Axial T2* FGRE
- Axial DWI b1000


Optional Series

- 3D Sag CUBE T2 FLAIR
- 3D Axial SWAN
- 3D Sag T1 MP-RAGE


2D Axial T2 FLAIR




3D Sag CUBE T2 FLAIR




Axial T2* FGRE




3D Axial SWAN



Axial DWI b1000



3D Sag T1 MP-RAGE



Take your imaging to the next level with SIGNA™,2



Use **AIR™ Recon DL** to reduce scan time by up to 50%



Use **AIRx™ Brain** to automatically place slices and reduce setup



These protocols were developed with the HNU; the **48ch Head Coil** will improve the available acceleration factors



Leverage integrated **cloud-based quantification** tools for white-matter, grey-matter and segmentation

1. Cogswell PM, et al. AJNR Am J Neuroradiol. 2022;43(9):E19-E35.

2. Some options are available for purchase
JB25532XX | May 2024

ARIA Monitoring Consensus Protocol - 3.0T

2D Axial T2
FLAIR

Axial T2*
FGRE

Axial DWI
b1000

3D Sag
CUBE T2
FLAIR

3D Axial
SWAN

3D Sag T1
MP-RAGE

Imaging Parameters

Imaging Mode	2D	2D	2D	3D	3D	3D
Pulse Sequence	FSE	Gradient Echo	Spin Echo EPI	CUBE T2 FLAIR	SWAN	MP-RAGE
Imaging Options	FC, TRF, Fast, ARC, T2flair	Fast, ARC	EPI, DIFF, Asset Classic	FAST, Z512, T2Prep, ARC, HS, IrP	FC, FAST, ARC, HS, Z512	Z512, ARC, HS, IrP, FAST
Phase Acceleration	2.0	2.0	2.0	2.0	2.0	2.0
Slice Acceleration	1.0	1.0	1.0	1.0	1.0	1.0
HyperSense Factor				1.2	1.2	1.2

Scan Range

FOV (cm)	22.0	22.0	24.0	22.0	20.0	25.6
Slice Thickness (mm)	4.0	4.0	4.0	0.6	1.0	0.5
Slice Spacing (mm)	0.4	0.4	0.4	-	-	-
Slice Resolution	-	-	-	50%	50%	50%

Acquisition Timing

Frequency	300	300	160	256	300	256
Phase	200	220	160	256	300	256
Phase FOV	0.9	0.9	1.0	1.0	0.9	1.0
NEX	1.0	1.0	1.0	1.0	0.7	1.0
Flow Comp Direction	Slice	-	-	-	-	-
Phase Image	-	-	-	-	Yes	-

Scan Timing

Flip Angle	Auto (160)	20.0	-	-	15.0	8.0
TE	130.0	15.0	Minimum	130.0	25.0	3.0
TR	9,000	800	7,500	6,000	Minimum	7.3
TI	Auto (2,474)	-	-	Auto (1,730)	-	1,000
# Echoes	1	1	1	1	4	1
Echo Train Length	23	-	-	230	-	-
Chem Sat	Fat	-	Fat	-	-	-
Receiver Bandwidth	62	50	250	50	31.25	31

The values in the table represent scan parameters developed with the 3.0T SIGNA™ Premier system using AIR™ Recon DL. Final parameters values may differ depending on your system's available software options.

3 Plane Localizer

3-Plane Localizer - 3-Plane Localizer	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	TE	80.0
	Pulse Sequence	Spin Echo	Number of Echoes	1
	Imaging Options	Seq, EDR, TRF, Fast, SS, ARC	TR	1000.0
	SCAN RANGE		Receiver Bandwidth	83.33
	FOV	24.0	IMAGE ENHANCE	
	Slice Thickness	10.00	Filter Choice	None
	Slice Spacing	5.0	USER CVS	
	Overlap Locations	0	User CV1	1.00
	ACQ TIMING		User CV Mask2	0
	Freq	256	MULTI-PHASE	
	Phase	192	Seperate Series	0
	Freq DIR	Unswap	Mask Phase	0
	# of Acq. Before Pause	0	Mask Pause	0
	Phase FOV	1.00	Preserve	0
	Auto Shim	Auto	DIFFUSION	
	Phase Correction	No	Recon All Images	On
	Excitation Mode	Selective	# Synthetic b-values	1
	FMRI		Synthetic b-value	1000.0;
	PSD Trigger	Internal	CONTRAST	
View Order	Bottom/Up	Contrast Yes/No	No	
# of Repetitions REST	0			
# of Repetitions ACTIVE	0			
SAT				
Tag Type	None			
TRICKS				
Pause On/Off	On			
Auto Subtract	0			
Auto SCIC	Off			

3-Plane Localizer - 3-Plane Localizer

2D Axial T2 FLAIR

Ax T2 FLAIR FS - Ax T2 FLAIR FS	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	Flip Angle	160.0
	Pulse Sequence	FSE-XL	TE	130.0
	Imaging Options	FC, EDR, TRF, Fast, ARC, T2flair	Number of Echoes	1
	Phase	2.00	TR	9000.0
	SCAN RANGE		TI	2473
	FOV	22.0	Echo Train Length	23
	Slice Thickness	4.00	Receiver Bandwidth	62.50
	Slice Spacing	0.4	IMAGE ENHANCE	
	Overlap Locations	0	Filter Choice	None
	Number of Slices	36	USER CVS	
	ACQ TIMING		User CV19	1.00
	Freq	300	User CV Mask2	128
	Phase	200	MULTI-PHASE	
	Freq DIR	A/P	Seperate Series	0
	Fat Shift DIR	Normal (A)	Mask Phase	0
	NEX	1.00	Mask Pause	0
	Phase FOV	0.90	Preserve	0
	Auto Shim	Auto	DIFFUSION	
	Phase Correction	No	Recon All Images	On
Flow Direction Compensation	Slice	Multi b-values	1000.0;	
RF Drive Mode	Preset	Multi NEX Values	1.0;	
Excitation Mode	Selective	# Synthetic b-values	1	
FMRI		Synthetic b-value	1000.0;	
PSD Trigger	Internal	CONTRAST		
View Order	Bottom/Up	Contrast Yes/No	No	
# of Repetitions REST	0			
# of Repetitions ACTIVE	0			
SAT				
Tag Type	None			
Fat/Water Saturation	Fat			
TRICKS				
Pause On/Off	On			
Auto Subtract	0			
Auto SCIC	Off			

Ax T2 FLAIR FS - Ax T2 FLAIR FS

Axial T2* FGRE

Ax T2* FGRE - Ax T2* FGRE	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	Flip Angle	20.0
	Pulse Sequence	Gradient Echo	TE	15.0
	Imaging Options	EDR, Fast, ARC	Number of Echoes	1
	Phase	2.00	TR	800.0
	SCAN RANGE		Receiver Bandwidth	50.00
	FOV	22.0	IMAGE ENHANCE	
	Slice Thickness	4.00	Filter Choice	None
	Slice Spacing	0.4	GATING/TRIGGER	
	Overlap Locations	0	Pause After Navigator	0
	Number of Slices	36	Prescan	
	ACQ TIMING		FMRI	
	Freq	300	PSD Trigger	Internal
	Phase	220	View Order	Bottom/Up
	Freq DIR	A/P	# of Repetitions REST	0
	NEX	1.00	# of Repetitions ACTIVE	0
	Phase FOV	0.90	SAT	
	Auto Shim	Auto	Tag Type	None
	Phase Correction	No	TRICKS	
	RF Drive Mode	Preset	Pause On/Off	On
Excitation Mode	Selective	Auto Subtract	0	
USER CVS		Auto SCIC	Off	
User CV Mask2	0	MULTI-PHASE		
MULTI-PHASE		Seperate Series	0	
Mask Phase	0	Mask Phase	0	
Mask Pause	0	Mask Pause	0	
Preserve	0	DIFFUSION		
DIFFUSION		Recon All Images	On	
Multi b-values	1000.0;	Multi b-values	1000.0;	
Multi NEX Values	1.0;	Multi NEX Values	1.0;	
# Synthetic b-values	1	# Synthetic b-values	1	
Synthetic b-value	1000.0;	Synthetic b-value	1000.0;	
CONTRAST		CONTRAST		
Contrast Yes/No	No	Contrast Yes/No No		

Ax T2* FGRE - Ax T2* FGRE

Axial DWI b1000

Ax DWI ALL b1000 - Ax DWI ALL b1000	IMAGING PARAMETERS		SCAN TIMING	
	Imaging Mode	2D	TE	Minimum
	Pulse Sequence	Spin Echo	Number of Echoes	1
	Imaging Options	EDR, Cla, EPI, DIFF, Asset	TR	7500.0
	Phase	2.00	Number of Shots	1
	SCAN RANGE		IMAGE ENHANCE	
	FOV	24.0	Filter Choice	None
	Slice Thickness	4.00	USER CVS	
	Slice Spacing	0.4	User CV0	1.00
	Overlap Locations	0	User CV5	1.00
	Number of Slices	36	User CV Mask2	0
	ACQ TIMING		MULTI-PHASE	
	Freq	160	Seperate Series	0
	Phase	160	Mask Phase	0
	Freq DIR	R/L	Mask Pause	0
	Phase FOV	1.00	Preserve	0
	Auto Shim	Auto	DIFFUSION	
	Phase Correction	Yes	Optimized TE	Yes
	RF Drive Mode	Preset	Diffusion Directions	All
	Excitation Mode	Selective	Number of Diffusion Directions	3
	FMRI		Number of T2 Images	1
	PSD Trigger	Internal	Dual Spin Echo	Off
	View Order	Bottom/Up	Recon All Images	On
	# of Repetitions REST	0	Multi b-values	1000.0;
# of Repetitions ACTIVE	0	Multi NEX Values	1.0;	
SAT		NEX For T2	1.00	
Tag Type	None	# Synthetic b-values	0	
Fat/Water Saturation	Fat	CONTRAST		
TRICKS		Contrast Yes/No	No	
Pause On/Off	On			
Auto Subtract	0			
Auto SCIC	Off			

Ax DWI ALL b1000 - Ax DWI ALL b1000

3D Sag CUBE T2 FLAIR

3D Sag T2 FLAIR Cube - 3D Sag T2 FLAIR Cube

3D Sag T2 FLAIR Cube - 3D Sag T2 FLAIR Cube

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	Cube T2 FLAIR
Imaging Options	EDR, Fast, ZIP512, T2P, ARC, HS, IrP
Phase	2.00
Slice	1.00
HyperSense	1.20

SCAN RANGE

FOV	22.0
Slice Thickness	0.60
Location per Slab	240
Overlap Locations	0
Number of Slices	1

ACQ TIMING

Freq	256
Phase	256
Freq DIR	S/I
Fat Shift DIR	Normal (S)
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Preset
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
Fat/Water Saturation	Fat

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

SCAN TIMING

TE	130.0
Number of Echoes	1
TR	6000.0
TI	1728
Echo Train Length	230
Receiver Bandwidth	50.00

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV5	1.00
User CV11	30.00
User CV34	100.00
User CV Mask2	10

MULTI-PHASE

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST

Contrast Yes/No	No
-----------------	----

3D Axial SWAN

3D Ax SWAN - 3D Ax SWAN

3D Ax SWAN - 3D Ax SWAN

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	SWAN
Imaging Options	FC, EDR, Fast, ARC, HS
Phase	2.00
Slice	1.00
HyperSense	1.20

SCAN RANGE

FOV	20.0
Slice Thickness	1.00
Location per Slab	160
Overlap Locations	0
Number of Slices	1

ACQ TIMING

Freq	300
Phase	300
Freq DIR	A/P
Phase FOV	0.90
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Preset
Excitation Mode	Selective
Phase Image	1

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
----------	------

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

SCAN TIMING

Flip Angle	15.0
TE	25.0
Number of Echoes	5
TR	Minimum
Receiver Bandwidth	31.25

IMAGE ENHANCE

Filter Choice	A
---------------	---

USER CVS

User CV22	1.00
User CV34	70.00
User CV Mask2	8

MULTI-PHASE

Seperate Series	0
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST

Contrast Yes/No	No
-----------------	----

3D Sag T1 MP-RAGE

3D Sag T1 MP-RAGE Fast HyperSense - 3D Sag T1 MP-RAGE Fast HyperSense

3D Sag T1 MP-RAGE Fast HyperSense - 3D Sag T1 MP-RAGE Fast HyperSense

IMAGING PARAMETERS

Imaging Mode	3D
Pulse Sequence	MP-RAGE
Imaging Options	EDR, Fast, ZIP512, ARC, HS, IrP
Phase	2.00
Slice	1.00
HyperSense	1.20

SCAN RANGE

FOV	25.6
Slice Thickness	0.50
Location per Slab	288
Overlap Locations	0
Number of Slices	1

ACQ TIMING

Freq	256
Phase	256
Freq DIR	S/I
Fat Shift DIR	Normal (S)
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	No
RF Drive Mode	Preset
Excitation Mode	Selective

FMRI

PSD Trigger	Internal
View Order	Bottom/Up
# of Repetitions REST	0
# of Repetitions ACTIVE	0

SAT

Tag Type	None
----------	------

TRICKS

Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

SCAN TIMING

Flip Angle	8.0
Number of Echoes	1
TI	1000
Receiver Bandwidth	31.25
Recovery Time	900

IMAGE ENHANCE

Filter Choice	None
---------------	------

USER CVS

User CV6	1.00
User CV Mask2	16384

MULTI-PHASE

Seperate Series	0
Trigger Delay without AV	0
Mask Phase	0
Mask Pause	0
Preserve	0

DIFFUSION

Recon All Images	On
# Synthetic b-values	1
Synthetic b-value	1000.0;

CONTRAST

Contrast Yes/No	No
-----------------	----